

Actigraphy measurement



- Relaxation-Activity-Profile for determination of circadian rhythm
- Recording of the motor activity for:
 - ▶ Sleep-Wake-Estimation
 - ▶ ADHS Diagnostics
 - ▶ Training / Sport / Rehabilitation
 - ▶ PLM recorder
- Tremor-Analysis (determination of frequency via FFT)

SOMNOwatch™ plus Acti


- ▶ Recording time up to 15 days
- ▶ Storage of raw data with 128 Hz
- ▶ Easy handling for patient and physician
- ▶ Patient marker for identification of relevant events during recording
- ▶ Built-in sensor for ambient light and body position



Sleep wake profile during 7 days

Recording of motor activity

Sleep/Wake estimation* and circadian rhythm



Sleep Wake Periods


Duration	Sleep Time	Wake Time	SPT	# Wake	Eff. (%)	Latency	WASO	MAI (h)	Σ Akt. (mg)
1 TIB: from 31.01 21:02:33 to 01.02 05:58:44									
08:56:10	08:05:44	00:50:26	08:49:14	4 (0,5)	90,6	00:06:56	00:43:30	5,1	16085
2 TIB: from 01.02 21:59:33 to 02.02 06:12:03									
08:12:30	07:48:03	00:24:27	08:09:03	3 (0,4)	95,0	00:03:27	00:21:00	4,7	4925
3 TIB: from 02.02 22:38:03 to 03.02 07:46:02									
09:09:59	07:30:30	01:39:29	07:40:30	3 (0,4)	81,9	00:03:27	01:36:02	5,7	51700
4 TIB: from 03.02 22:54:33 to 04.02 07:52:03									
08:57:30	06:57:30	02:00:00	07:30:00	6 (0,9)	77,7	00:04:57	01:55:03	5,2	58050
5 TIB: from 04.02 22:45:33 to 05.02 09:23:33									
10:38:00	10:25:33	00:12:27	10:28:03	1 (0,1)	98,0	00:09:57	00:02:30	5,7	6380
6 TIB: from 05.02 23:12:03 to 06.02 09:35:32									
10:23:29	09:08:00	01:15:29	10:11:30	5 (0,5)	87,9	00:02:57	01:12:32	4,8	14964
7 TIB: from 06.02 22:05:03 to 07.02 07:49:03									
09:44:00	07:42:30	02:01:30	08:21:00	5 (0,6)	79,2	00:02:27	01:59:03	6,2	58400
Average - TIB									
09:25:57	08:13:58	01:11:58	08:44:11	3 (0,1)	87,2	00:04:52	01:07:05	5,3	30072

- For detection of Sleep/Wake-rhythm: Application to the non-dominant arm
- Recording of activity in epochs of 1- 120s and display in a actigraphy profile
- Determination of the relevant recording time and day/night with the built-in light sensor

ADHS diagnostics



inactive [min]	Acti [mg]	active [min]	inactive [min]
School von 30.05 08:44:59 bis 30.05 14:44:59	182	186,6 (51,8%)	173,4 (48,2%)
Leisure von 30.05 14:44:59 bis 30.05 22:19:59	149	200,0 (43,9%)	255,0 (56,1%)
Sleep von 30.05 22:20:00 bis 31.05 08:44:59	42	35,1 (5,6%)	589,9 (94,4%)
School von 31.05 08:44:59 bis 31.05 14:45:03	207	190,1 (52,8%)	170,0 (47,2%)
Leisure von 31.05 14:45:04 bis 31.05 22:20:00	158	224,1 (49,3%)	230,8 (50,7%)
Sleep von 31.05 22:20:00 bis 01.06 08:44:59	53	45,2 (7,2%)	579,8 (92,8%)

School 12:00:04	196 194 ± 35	188,4 188,3 ± 4,9	171,7 171,7 ± 4,9
Leisure 15:09:56	153 154 ± 13	212 212 ± 34,1	
Sleep 20:49:59	48 48 ± 16	40,2 40,2 ± 14,4	
Total	5651,7		



- For ADHS diagnostics pre-selected periods of time can be allocated to specific activities, such as school, breaks and recreation

Tremor Analysis



- Determination of low- and high frequency tremor and amplitude
- Graphic presentation of tremor frequency and intensity
- Calculation of average tremor frequency and duration
- Detailed display of frequency range (adjustable limits)

7 internal channels

Body position, 3 activity sensors (x, y, z-axis, magnitude), ambient light, patient marker
Data recording adjustable sampling rate up to 128/s, small file sizes thanks to compressed data storage
Power supply Li-Ion battery, 630mAh, rechargeable
Size and weight 45 mm diameter x 16 mm (incl. battery), 30 g
Device functions Recording duration up to 15 days Battery charging and data transfer via docking station programmable start and endtimes, delayed auto start splashproof

* Dick, R., et al., AASM standards of practice compliant validation of actigraphic sleep analysis from SOMNOwatch versus polysomnographic sleep diagnostics shows high conformity also among subjects with sleep disordered breathing. Physiological measurement, 2010. 31(12): p. 1623-33.