


Actigraphy measurement



- Relaxation-Activity-Profile for determination of circadian rhythm
- Recording of the motor activity for:
 - ▶ Sleep-Wake-Estimation
 - ▶ ADHS Diagnostics
 - ▶ Training / Sport / Rehabilitation
 - ▶ PLM recorder
- Tremor-Analysis (determination of frequency via FFT)

SOMNOwatch™ plus Acti


- ▶ Recording time up to 15 days
- ▶ Storage of raw data with 128 Hz
- ▶ Easy handling for patient and physician
- ▶ Patient marker for identification of relevant events during recording
- ▶ Built-in sensor for ambient light and body position



Sleep wake profile during 7 days

Recording of motor activity

Sleep/Wake estimation* and circadian rhythm




Sleep Wake Periods

Duration	Sleep Time	Wake Time	SPT	# Wake	Eff. (%)	Latency	WASO	MAI (h)	Σ Akt. (mg)
1 TIB: from 31.01.21:02:33 to 01.02.06:58:44									
08:56:10	08:05:44	00:50:20	08:49:14	4 (0,5)	99,8	00:06:58	00:43:30	5,1	16085
2 TIB: from 01.02.21:59:33 to 02.02.06:12:03									
08:12:30	07:48:03	00:24:27	08:09:03	3 (0,4)	95,0	00:03:27	00:21:00	4,7	4925
3 TIB: from 02.02.22:38:03 to 03.02.07:48:02									
09:09:59	07:39:30	01:29:29	07:40:30	3 (0,4)	81,9	00:03:27	01:36:02	5,7	51700
4 TIB: from 03.02.22:54:33 to 04.02.07:02:03									
08:57:30	08:57:30	02:00:00	07:30:00	6 (0,9)	77,7	00:04:57	01:55:03	5,2	58050
5 TIB: from 04.02.22:48:33 to 06.02.09:23:33									
10:38:00	10:25:53	00:12:27	10:28:03	1 (0,1)	98,0	00:09:57	00:02:30	5,7	8380
6 TIB: from 05.02.23:12:03 to 06.02.09:35:32									
10:23:29	09:08:00	01:15:29	10:11:30	5 (0,5)	87,9	00:02:57	01:12:32	4,8	14964
7 TIB: from 06.02.22:05:03 to 07.02.07:45:03									
09:44:00	07:42:30	02:01:30	08:21:00	5 (0,8)	79,2	00:02:27	01:59:03	6,2	58400
Average - TIB									
09:25:57	08:13:58	01:11:58	08:44:11	3 (0,1)	87,2	00:04:52	01:07:05	5,3	30072

- For detection of Sleep/Wake-rhythm: Application to the non-dominant arm
- Recording of activity in epochs of 1- 120s and display in a actigraphy profile
- Determination of the relevant recording time and day/night with the built-in light sensor



ADHS diagnostics

inactive [min]	Acti [mg]	active [min]	inactive [min]
School von 30.05.08:44:50 bis 30.05.14:44:50			
180	186,6 (51,8%)	173,4 (48,2%)	
Leisure von 30.05.14:44:50 bis 30.05.22:19:59			
149	200,0 (43,9%)	255,0 (56,1%)	
Sleep von 30.05.22:20:00 bis 31.05.08:44:59			
42	35,1 (5,6%)	589,9 (94,4%)	
School von 31.05.08:44:59 bis 31.05.14:45:00			
207	190,1 (52,8%)	170,0 (47,2%)	
Leisure von 31.05.14:45:04 bis 31.05.22:20:00			
158	224,1 (49,3%)	230,8 (50,7%)	
Sleep von 31.05.22:20:00 bis 01.06.08:44:59			
53	45,2 (7,2%)	579,8 (92,8%)	



- For ADHS diagnostics pre-selected periods of time can be allocated to specific activities, such as school, breaks and recreation


Tremor Analysis



- Determination of low- and high frequency tremor and amplitude
- Graphic presentation of tremor frequency and intensity
- Calculation of average tremor frequency and duration
- Detailed display of frequency range (adjustable limits)

7 internal channels
Body position, 3 activity sensors (x, y, z-axis, magnitude), ambient light, patient marker
Data recording
adjustable sampling rate up to 128/s, small file sizes thanks to compressed data storage
Power supply
Li-Ion battery, 630mAh, rechargeable
Size and weight
45 mm diameter x 16 mm (incl. battery), 30 g
Device functions
Recording duration up to 15 days
Battery charging and data transfer via docking station
programmable start and endtimes, delayed auto start
splashproof

* Dick, R., et al., AASM standards of practice compliant validation of actigraphic sleep analysis from SOMNOwatch versus polysomnographic sleep diagnostics shows high conformity also among subjects with sleep disordered breathing. Physiological measurement, 2010. 31(12): p. 1623-33.



SOMNOmedics GmbH · Am Sonnenstuhl 63 · 97236 Randersacker · Germany
Tel: +49-931 / 359094-0 · info@somnomedics.eu · www.somnomedics.eu