


SOMNOmedics

Respiratory Screening Designed for maximum comfort

SOMNOwatch™ plus RESP

- 
- Breathing / Flow (nasal cannula)
 - Effort Thorax & Abdomen (RIP)
 - Snoring
 - SpO₂
 - Pulse rate
 - Body position
 - Movement (Sleep/Wake)
 - CPAP/BiPAP-Pressure
 - Ambient Light (TIB)
 - Patient Marker

- + 80hrs recording time by data compression
- + Programmable alarm sound
- + Built in rechargeable Li-Ion battery
- + Fast data transfer by USB Port
- + Recorded data directly sent via E-Mail
- + Calculation of AHI & RDI
- + Low maintenance costs