

SOMNOwatch™ PLUS

NEED YOUR WEARABLE TO WORK HARDER? OURS DOES.

Research Grade Sensors

With programable software allowing extensive customisation

7 Internal Channels

Body Position (Acc for 3 axis), Ambient light, patient marker with acoustic tone

High sampling rates

256/32 Hz

Data Compression

Almost no limit to time duration

External Sensor Compatible

Up to 8 channels inc, BP, Resp, 6 channel EEG, long term ear temp and much much more

Tremor Analysis

Intentional Tremor < 5Hz in the FFT of the x-axis

