

Gold Cup Electrodes

The following product codes can be used to order the Gold Cup Electrodes




Product	Item Nr.
Gold Cup Electrodes, 75 cm (30") wire, 10 pcs.	SEN710
Gold Cup Electrodes, 100 cm (40") wire, 10 pcs.	SEN711
Gold Cup Electrodes, 150 cm (60") wire, 10 pcs	SEN712
Gold Cup Electrodes, 200 cm (80") wire, 10 pcs.	SEN713
Gold Cup Electrodes, 250 cm (98") wire, 10 pcs.	SEN714



Order easily and conveniently:

CONTACT

- phone +49 (931) 359 094-0
(Mo.-Fr.: 08:00-17:00)
- fax +49 (931) 359 094-49
- e-Mail sales@somnomedics.de
- online www.somnomedics.eu



© 2018 SOMNOmedics GmbH, Germany



Subject to change (without notice)

9_5_01-2019_EN_Rev_0

Gold Cup Electrodes

- ▶ For a constant and reliable EEG recording
- ▶ Easier handling
- ▶ Exceptionally robust
- ▶ Outstanding Price / Quality ratio

Can be used for EEG, PSG, EMG, EP, OEG,
ECG and LTM (Long Term Monitoring).



Gold Cup Electrodes

The special shape of these Gold Cup Electrodes allows easy placement of the electrodes on the scalp (EEG) or on the patients face (EOG).

They also provide a constant and reliable signal during EEG, PSG, EMG, ECG, EOG and Video-EEG-Long-Term-Monitoring (LTM). Our low impedance electrodes are equipped with a special low-noise cable.

Properties:

- ◆ Excellent recording quality
- ◆ Highly flexible wire for easy application, cleaning and disinfection
- ◆ Ultra low noise
- ◆ Bend protection to ensure durability
- ◆ Various cable lengths
- ◆ Compatible with various EEG and PSG systems
- ◆ Can be used with most popular gels, pastes and adhesives
- ◆ Color coded lead wires for easy handling (10 colors)
- ◆ 10 mm diameter
- ◆ 1.5 mm female DIN42802 connector
- ◆ 10 electrodes per package

Compatibility

of Gold Cup Electrodes:

Gold Cup Electrodes are compatible with all medical devices with an 1.5 mm input socket including:

- | | |
|-----------------------|----------------|
| ◆ SOMNOmedics | ◆ Grass |
| ◆ Embla | ◆ Xltek |
| ◆ Compumedics | ◆ Lifelines |
| ◆ Natus | ◆ Nicolet |
| ◆ Cadwell | ◆ NicoletONE |
| ◆ Dr. Fenyves und Gut | ◆ Nihon Kohden |
| ◆ Teca Nervus | ◆ Weinmann |
| ◆ Philips Respironics | ◆ Nox |
| ◆ Löwenstein Medical | ◆ Bio-logic |



Gold Cup Electrodes

lic2® Electrode Cream
Pcs: 10 Tubes
Item No.: SEN035