

# SOMNOtouch™ PSG

## OUR BIGGEST IDEAS. NOW IN OUR SMALLEST PSG DEVICE

### Headbox with Continuous Impedance Check

Lower noise, higher signal quality and easy confirmation of proper electrode hook-up

### AASM Compliant Recording (up to 25 Channels)

Headbox with EEG, EOG, EMG. Combi sensor for PLM, ECG & effort. 11 Internal channels inc: SpO<sub>2</sub>, pulse rate & plethysmogram

### Real-time Telemetric Data Transfer

On a PC, Smartphone or Tablet

ACTUAL  
SIZE

### Cuffless Continuous Blood Pressure Measurement

A new paradigm in measuring blood pressure in a sleep lab – without disturbing sleep (optional upgrade)

### Intelligent Connect

Automatic sensor recognition – making montages a thing of the past

### Small Compact Robust Design

Credit card footprint means unparalleled comfort. Touch Screen display for PC independent set-up & signal quality check at bedside




The **SOMNOtouch™ PSG** combines the comfort, function and power of much larger sleep diagnostic devices in a small perfectly formed size. At only 64 g and 84 mm × 55 mm it is one of the lightest yet powerful touch screen ambulatory/stationary polysomnography devices on the market today.

## AUTOMATIC SOFTWARE & FIRMWARE UPDATES

### Colour Touch Screen

- ◆ Start the data recording on the device
- ◆ Signal check on the screen
- ◆ Zoom to a selectable time base
- ◆ Easy to use interface
- ◆ Clear icons to easily identify attached sensors



### Versatility – In a compact design

- ◆ Respiratory screening, plus recording of EEG, EMG & EOG signals from the headbox and of the ECG & PLM signals with a combi-sensor
- ◆ Minimal costs per measurement
- ◆ Configurable across reimbursement systems



PTZ Camera incl. Infrared (optional), Item no: VID150A

### Intelligent Connect

- ◆ Automatic sensor recognition, no montages required
- ◆ Clear icons to easily identify functions/sensors
- ◆ AASM compliant headbox with continuous impedance recording

Never miss a new feature or update

### Video & Accessories (optional)

- ◆ Now with on-line recording capabilities
- ◆ Powered by our powerful DOMINO software
- ◆ Complete range of stationary video solutions
- ◆ Comfortable Bluetooth receiver (BTR062)
- ◆ Include external devices via our ADC board

### Docking Station

- ◆ Fast recharging for PSG recording
- ◆ Fast data transfer (Approx. 90s for a PSG)
- ◆ Virtual dockingstation to display device status on PC
- ◆ Available in a 3 port design for bigger installations

